NADAL® pH Test

REF 194002

DE Gebrauchsanweisung	2	PT Instruções de Utilização	20
EN Instructions for use	5	CZ Návod k použití	23
FR Instructions d'utilisation	8	FI Käyttöohje	26
ES Instrucciones de uso	11	SE Användarinstruktioner	29
IT Istruzioni per l'uso	14	Symbols	35
PL Sposób użycia	17	Our Teams	36

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14. Referenzen

1. Boskey E.R. et al.: Acid Production by Vaginal Flora In Vitro is Consistent with the Rate and Extent of Vaginal Acidification, Infect Immun (1999)
2. Boskey E.R. et al.: Origins of vaginal acidity: high D/L lactate ratio is consistent with bacteria being the primary source, Human Reproduction (2001)
3. Kiss et al., Prospective randomized controlled trial of an infection screening programme to reduce the rate of preterm delivery, BMJ (2004)
4. Ferris D. G. et al.: Variability of vaginal pH determination by patients and clinicians, JABFM (2006)
5. Heinze T. et al.: Determination of vaginal pH by pH indicator strip and by pH micro electrode, J Perinat Med (1989)

Rev. 0, 2014-11-26 OM/UV

1. Intended Use

The NADAL® pH Test is for measuring the pH value of vaginal fluid. The vaginal pH value is a useful screening parameter for disturbances of the vaginal microflora. Abnormal pH values are frequently an indicator for clinical or sub-clinical vaginal infections.


The intended uses of the NADAL® pH Test are as follows:

1. The test can be used as an aid in the diagnosis of vaginal infections or in the identification of conditions that are supportive for the manifestation of vaginal infections in pregnant and non-pregnant premenopausal women.
2. The test can be used for regular screening of the vaginal pH in pregnant women in order to increase the chances of early detection and treatment of vaginal infections, thus to reduce the risk of premature birth.

The test is designed for *in-vitro* diagnostic use by professionals. It should not be used without proper supervision and is not intended for over the counter sale to the general public.

2. Introduction

A physiological acidic vaginal pH value of 3.8 to 4.4 protects the vagina from colonization by pathogens thus decreasing the incidence of vaginal infections. Mainly Lactobacilli as part of the natural vaginal microflora contribute to the acidic pH value by the degradation of glycogen released by cells of the upper vaginal wall to lactic acid under the anaerobic conditions present in the mucosa of the vagina.

An increase of the vaginal pH value can be caused by an untypical bacterial colonization of the vagina (bacterial vaginosis), vaginal infections or STDs. But also other conditions like PROM (Premature rupture of membranes) events, estrogen deficiency, postoperative infected wounds, excessive intimate care or antibiotic therapy might affect the pH value.

The pH test can be used as an aid in diagnosing vaginal infections in women showing typical symptoms of infection like inflammation, itching, increased discharge etc. Regular

testing of vaginal pH value even in non-symptomatic women will help to identify conditions that support the proliferation of bacteria that might lead to gynecological complaints and complications including sterility.

For pregnant women a major benefit of regular vaginal pH testing is in the reduction of the risk of premature birth. It is known that some pathogens contribute to early contractions or might migrate upwards the vagina causing severe infections that might lead to complications in pregnancy and birth. Regular pH screening will contribute to the identification of these conditions allowing an early treatment of the patients.

3. Test Principle

The NADAL® pH Test is a reliable, hygienic method for determining the vaginal pH value.

It consists of a vaginal applicator of approx. 8 cm in length with an integrated pH measurement strip which contains a pH indicator with an optimal range of resolution in the spectrum of interest (MERCCK pH indicator rods, pH 4.0-7.0, not bleeding, 1.09542.0001). The applicator consists of a round handle area and a cone-shaped insertion stick of approx. 5 cm in length. On one side of the tip of the insertion stick there is a window-like opening, where the pH measurement strip is located (pH measurement area).

The round handle allows safe handling of the applicator. The test is inserted approx. 1-3 cm into the vagina and the pH measurement area is gently pressed against the back of the vaginal wall. As soon as the pH measurement area gets moistened with vaginal fluid the colour change occurs. The vaginal applicator is removed from the vagina and the pH value is read by assigning the colour to a value on the provided colour scale.

4. Reagents and Materials Supplied

- 50 hygienically packaged NADAL® pH Test applicators for the measurement of vaginal pH
- 1 package insert
- 1 pH colour scale for reading the result

5. Additional Materials Required

No additional materials are required.

6. Storage & Stability

Test kits should be stored between 2-30°C. Tests should remain in the sealed foil packaging until ready to use. Do not freeze. Do not expose tests and colour scale to direct sunlight. Under these conditions tests remain stable until the expiration date printed on the pouch.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Applicators for single use only.
- Do not use tests after the stated expiry date.
- Do not use tests if pouch is damaged.
- Do not touch the pH measurement area of the applicator to prevent contamination.
- Sample material is potentially infectious. Standard guidelines for handling infectious agents and chemical reagents should be observed throughout all procedures. Protective garments (gloves, laboratory coats, eye protection) are recommended. Material that has been in

contact with potentially infectious sample material should be disposed of according to local regulations.

- Observe instructions for specimen collection in order to avoid conditions that might temporarily affect vaginal pH value leading to false results.
- Observe instructions for result interpretation. Computer simulations suggest that colour blindness should not affect correct reading but we recommend supervision of colour blind interpreters until correct reading ability is proven. Results should be read immediately while test is still moist.

8. Specimen Collection and Storage



The following conditions are known to affect the vaginal pH temporarily and should be avoided before the pH measurement:

- Remaining urine at the vaginal entrance – it is recommended to wait 15 minutes after urination and to ensure by wiping to remove remaining urine before performing the test.
- Seminal fluid - measurement should take place no earlier than 12 hours after sexual intercourse.
- Vaginal medical products (vaginal suppositories, creams, gels, etc.) – measurement should be done no earlier than 12 hours after last use.
- Menstruation – measurement should not be done during menstruation. It is recommended to wait 3-4 days after menstruation before using the test.
- Intense intimate care: Genital area should not have been washed immediately before taking the measurement. Soaps, shower gels, etc. might affect the pH. Also vaginal douches are not recommended.

9. Test Procedure

If tests have been stored refrigerated, allow them to come to room temperature before starting the test procedure.

1. Remove the pH test from the sealed pouch, holding it at the round handle. Do not touch the pH measurement zone to avoid contamination. After having been removed from the protective pouch the test should be used immediately.
2. Insert the applicator into the vagina in a way that the pH measurement area comes into contact with the back wall of the vagina. Ensure to minimize contact with other liquids during this process. Insertion depth should be not deeper than 2-3 cm (approx. 1 inch) and should not reach beyond the round handle.
3. Leave the applicator in the vagina for approximately 10 seconds applying slight pressure to the vaginal wall to allow moistening of the pH measurement area.
4. Withdraw the applicator from the vagina and immediately start the interpretation of the result

10. Result Interpretation

Read the result immediately while the pH measurement area is still moist.

1. Compare the colour in the pH measurement area with the provided colour scale for the semiquantitative reading.
2. The colour that matches the colour of the pH measurement area best, represents the pH level.

Important: If the pH measurement area does not display a homogeneous colour, it might be that the area has not been sufficiently moistened with vaginal fluid. If this is the case, the measurement can be repeated using the same vaginal applicator. Alternatively, you can use a new test.

pH value in normal range

Vaginal pH levels in the range of 4.0 to 4.4 are within the normal/healthy range. The colour of the indicator is yellowish.

pH value in atypical range

Vaginal pH levels over 4.4 are outside of the normal/healthy range. At these levels the colour of the indicator is darker and goes into a brown-green-blue shade.

As some conditions might temporarily affect the vaginal pH, it is recommended to repeat testing a few hours later or the next morning. If the pH value remains elevated, further investigations/medical actions should be taken into consideration.

11. Quality Control

External controls are not provided with this kit. It is recommended to verify test performance as required by internal lab procedures. Solutions of known pH value might be used as controls.

12. Limitations

- The applicator determines only the vaginal pH value and not the presence of infection. As with all diagnostic tests, a definite clinical diagnosis should not be based on the results of the test alone, but should only be made by the physician after all clinical and laboratory findings have been evaluated.
- In menopausal/postmenopausal women the vaginal microflora changes due to low estradiol levels. This frequently results in elevated vaginal pH levels.
- The pH indicator used in the pH measurement area can only resolve the pH range printed on the colour scale. Higher and lower pH values will be shown as the highest and the lowest value on the scale.
- It should be kept in mind that there are conditions that influence the pH value temporarily (see Sample Collection and Interference) and to ensure that patients have avoided such conditions.

13. Performance Characteristics

Measurement Range and analytical sensitivity

The measurement range of the NADAL® pH Test is between 4.0 and 7.0. Best resolution is found for the clinically important range between pH 4.4 and 5.5 allowing the user to discriminate between normal and atypical vaginal pH values. Lower values than pH 4.0 will be shown with the value of 4.0. Higher values than pH 7.0 will be shown with the value of 7.0.

Accuracy

Solutions of known pH were tested in 10 fold determinations with the NADAL® pH Tests in order to see if results could be aligned to the corresponding value of the colour scale. For the critical pH of 4.4 the standard deviation was calculated to be $\text{pH } 4.4 \pm 0.11$. So the test is capable of discriminating samples with normal and with atypical pH values with sufficient accuracy.

Clinical performance

The suitability of the pH indicator for measuring vaginal pH has been proven in a publication by Heinz et al (1989) by comparison against pH measuring with a microelectrode placed directly in the vagina.

Results from blind testing of solutions of known pH value (pH electrode) were evaluated regarding their physiological evidence. Specimens up to pH 4.4 were classified as positive samples and described as normal vaginal environment. Above pH 4.4 samples were classified as negative. These samples described an atypical vaginal environment with an increased risk of infection.

		pH electrode		
		+	-	Total
NADAL® pH Test	+	54	20	74
	-	0	127	127
	Total	54	147	201

Sensitivity: > 99.9%

Specificity: 86.4%

Positive predictive value (PPV): 73%

Negative predictive value (NPV): > 99.9 %

Accuracy: 90 %

Interference Testing

Acetaminophen, Ascorbic Acid, Salicylic Acid, Atropine, Caffeine, Glucose, Hemoglobin, Ampicillin, EDTA, Ephedrine, Ibuprofen, Ethanol, Phenylpropylamine, Tetracycline, Human Albumine, Bilirubin, Creatinine, Estradiol, Estriol, 3-Hydroxybutyric Acid, Pregnanediol, Progesterone, Riboflavin, Uric Acid, Urea, Lactic acid, Glycerol, Acetic acid were tested for their interference with result generation. As some of the substances influenced the pH, pH values were adjusted after the addition of the substances to the desired values. Except for hemoglobin none of the substances interfered with the result formation at pH values of 4.0; 4.4; 4.7; 5.0; 5.3; 5.5; 5.8 and 7.0.

For hemoglobin containing solutions the colour of hemoglobin prevented a correct reading at visible concentrations.

14. References

1. Boskey E.R. et al.: Acid Production by Vaginal Flora In Vitro is Consistent with the Rate and Extent of Vaginal Acidification, *Infect Immun* (1999)
2. Boskey E.R. et al.: Origins of vaginal acidity: high D/L lactate ratio is consistent with bacteria being the primary source, *Human Reproduction* (2001)
3. Kiss et al., Prospective randomized controlled trial of an infection screening programme to reduce the rate of preterm delivery, *BMJ* (2004)
4. Ferris D. G. et al.: Variability of vaginal pH determination by patients and clinicians , *JABFM* (2006)
5. Heinze T. et al.: Determination of vaginal pH by pH indicator strip and by pH micro electrode, *J Perinat Med* (1989)

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14. Referencias

1. Boskey E.R. et al.: Acid Production by Vaginal Flora In Vitro is Consistent with the Rate and Extent of Vaginal Acidification, Infect Immun (1999)
2. Boskey E.R. et al.: Origins of vaginal acidity: high D/L lactate ratio is consistent with bacteria being the primary source, Human Reproduction (2001)
3. Kiss et al., Prospective randomized controlled trial of an infection screening programme to reduce the rate of preterm delivery, BMJ (2004)
4. Ferris D. G. et al.: Variability of vaginal pH determination by patients and clinicians , JABFM (2006)
5. Heinz T. et al.: Determination of vaginal pH by pH indicator strip and by pH micro electrode, J Perinat Med (1989)

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14. Bibliografia

- Boskey E.R. et al.: Acid Production by Vaginal Flora In Vitro is Consistent with the Rate and Extent of Vaginal Acidification, Infect Immun (1999)
- Boskey E.R. et al.: Origins of vaginal acidity: high D/L lactate ratio is consistent with bacteria being the primary source, Human Reproduction (2001)
- Kiss et al., Prospective randomized controlled trial of an infection screening programme to reduce the rate of preterm delivery, BMJ (2004)
- Ferris D. G. et al.: Variability of vaginal pH determination by patients and clinicians, JABFM (2006)
- Heinze T. et al.: Determination of vaginal pH by pH indicator strip and by pH micro electrode, J Perinat Med (1989)

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14. Referências

1. Boskey E.R. et al.: Acid Production by Vaginal Flora In Vitro is Consistent with the Rate and Extent of Vaginal Acidification, Infect Immun (1999)
2. Boskey E.R. et al.: Origins of vaginal acidity: high D/L lactate ratio is consistent with bacteria being the primary source, Human Reproduction (2001)
3. Kiss et al., Prospective randomized controlled trial of an infection screening programme to reduce the rate of preterm delivery, BMJ (2004)
4. Ferris D. G. et al.: Variability of vaginal pH determination by patients and clinicians, JABFM (2006)
5. Heinze T. et al.: Determination of vaginal pH by pH indicador strip and by pH micro electrode, J Perinat Med (1989)

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14. Referenser

1. Boskey E.R. et al.: Acid Production by Vaginal Flora In Vitro is Consistent with the Rate and Extent of Vaginal Acidification, Infect Immun (1999)
2. Boskey E.R. et al.: Origins of vaginal acidity: high D/L lactate ratio is consistent with bacteria being the primary source, Human Reproduction (2001)
3. Kiss et al., Prospective randomized controlled trial of an infection screening programme to reduce the rate of preterm delivery, BMJ (2004)
4. Ferris D. G. et al.: Variability of vaginal pH determination by patients and clinicians , JABFM (2006)
5. Heinze T. et al.: Determination of vaginal pH by pH indicator strip and by pH micro electrode, J Perinat Med (1989)

Rev. 0, 2014-11-26 JB

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conformité aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consultense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	In-vitro-Diagnostika	In-vitro diagnostic medical device	Dispositif médical de diagnostic in vitro	Producto sanitario para diagnóstico in vitro	Dispositivo medico-diagnostico in vitro	Tylko do diagnostyki in vitro
	Temperaturbegrenzung	Temperature limitation	Limites de température	Límite de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Code du lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour "n" tests	Suficiente para <n> utilizaciones	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń

Symbol	Português	Český	Suomi	Svenskt	Nederlands	Dansk
	Conformidade com as normas europeias	CE certifikát	CE-merkitty	CE-märkning	CE-markering	CE-mærkning
	Consultar as instruções de utilização	Viz návod k použití	Katso käyttöohjetta	Läs bruksanvisningen	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen
	Dispositivo médico para diagnóstico in vitro	Diagnostický zdravotnický prostředek in vitro	In vitro - diagnostiikkaan tarkoitettu lääkinällinen laite	Medicinteknisk produkt avsedd för in vitro-diagnostik	Medisch hulpmiddel voor in-vitrodiagnostiek	Medicinsk udstyr til in vitro-diagnostik
	Limites de temperatura	Teplotní omezení	Lämpötilarajat	Temperaturbegränsning	Temperatuurlimiet	Temperaturbegränsning
	Código do lote	Kód šarže	Eräkoodi	Satsnummer	Code van de partij	Batchkode
	Não reutilizar	Pro jednorázové použití	Kertakäyttöinen	Får inte återanvändas	Niet opnieuw gebruiken	Må ikke genbruges
	Prazo de validade	Spotřebujte do	Käytettävä viimeistään	Används före	Houdbaar tot	Udløbsdato
	Número de catálogo	Katalogov číslo	Luettelonumero	Listnummer	Catalogus nummer	Bestil lingsnummer
	Fabricante	Výrobce	Valmistaja	Tilverkare	Fabrikant	Fabrikant
	Suficiente para <n> test	Dostačuje pro <n> testů	Lukumäärä <n> test	Räcker till <n> test	Volgende voor <n> test	Tilstrækkeligt til <n> test

Our Teams

Germany:

Regensburg

Tel: +49 941 290 10-0
Fax: +49 941 290 10-50

Moers

Tel: +49 2841 99820-0
Fax: +49 2841 99820-1

Austria:

Tel: +49 941 290 10-29
Free Tel: 0800 291 565
Fax: +49 290 10-50
Free Fax: 0800 298 197

UK & Ireland:

Tel: +49 941 290 10-18
Free Tel – UK: 0800 234 1237
Free Tel – IRE: 1800 555 080
Fax: +49 290 10-50

France:

France Tel: 0800 915 240
France Fax: 0800 909 493

Switzerland

Swiss Tel: 0800 564 720
Swiss Fax: 0800 837 476

Belgium

Belgium Tel: 0800 718 82
Belgium Fax: 0800 747 07

Luxembourg

Lux. Tel: 800 211 16
Lux. Fax: 800 261 79

Spain:

Tel: +49 941 290 10-759
Free Tel: 900 938 315
Fax: +49 941 290 10-50
Free Fax: 900 984 992

Italy:

Tel: +49 941 290 10-34
Fax: +49 941 290 10-50

Poland:

Tel: +49 941 290 10-44
Free Tel: 00 800 491 15 95
Fax: +49 941 290 10-50
Free Fax: 00 800 491 15 94

Portugal:

Tel: +49 941 290 10-735
Tel. Verde: 800 849 230
Fax: +49 941 290 10-50
Fax Verde: 800 849 229

Netherlands:

Tel: +31 30 75 600
Free Tel: 0800 0222 890
Fax: +31 70 30 30 775
Free Fax: 0800 024 9519

Nordic countries (Finland, Norway, Sweden, Denmark):

Tel: +31 703075 607
Free Tel: +45 80 88 87 53
Tax: +31 703030 775

Laboratory Diagnostics Team:

Tel: +49 941 290 10-40
Fax: +49 941 290 10-50

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany

www.nal-vonminden.com • info@nal-vonminden.com

Fon: +49 2841 99820-0 • Fax: +49 2841 99820-1

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nal von minden GmbH