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- 3. Lecture des Prick-tests**
 - Mark on the test sheet the name (Surname and first name) of the patient, add the doctor's seal or stamp, add the date and the time of the prick tests.
 - The reading** is best made, on average, after 20 minutes.
 - Clean the skin** thoroughly with alcohol, taking care not to mix the solutions by sweeping, to remove any mark of the glycerol saline diluents which could impede marking with the felt tip pen.
 - Use a **very fine felt tip pen** to draw an outline of the weal. The presence of the papule can be estimated by touch and the outline can be viewed precisely by use of an oblique light.
 - The test report** is gathered by laying a transparent self-adhesive tape on the skin sketching the outline of the weal down on the test sheets according to the test point.
- 4. Interpretation of results**
 - It is possible when the **skin reactivity** meets two conditions :
 - Positive control response: The weal diameter of a prick test with 9% codeine phosphate (or 10 mg/ml histamine hydrochloride) must be larger than or equal to 3 mm.
 - Negative control response: The weal diameter of a prick test with the glycerol saline diluent must be less than or equal to 3 mm.
 - The test is **positive** when :
 - The weal diameter is higher than 3 mm or 3 mm larger than the negative control one
 - The test is negative when the weal response is less than 3 mm.

Adverse reactions:
Risk of light bleeding at the site. In order to limit this, reduce the applied pressure.
Any serious incident that has occurred in relation to Stallerpoint should be reported to Stallergenes (service.pharmacovigilance@stallergenesreer.com) and the competent authority of the Member State in which the user and/or patient is established.

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Edition : 02/2022

Form and presentation:
Disposable Sterile Polymethacrylate Micropuncture Point. Stallerpoint is a medical device.

Box contents:
Large size: a box containing 1000 Stallerpoint® - 100 thin slabs of 10 Stallerpoint®.
Small size: a box containing 100 Stallerpoint® - 10 thin slabs of 10 Stallerpoint®.

Name and address of the manufacturer
STALLERGENES
6, rue Alexis de Tocqueville
F-92160 Antony
FRANCE
Tél.: (33) 1 55 59 20 00

Stallerpoint™ is used in association with a Prick-test allergen extract solution as part of the diagnostic test for immediate type hypersensitivity on the skin and for the standardization of allergen extracts. The device is used for children, adolescents and adults.

Methodology of the Prick-tests:
Stallerpoint is intended to be used by a healthcare professional trained to skin prick testing.

- Preparation of the Prick-tests**
 - Prepare a sheet for the collection of results.
 - Set out the vial with the Prick-tests on the display unit in the same order as on the test sheets.
 - Skin preparation:** Dry clean the forearm skin or the back with an alcohol-soaked cotton or alcohol wipe and let dry.
 - Mark the spots** for each of the Prick tests to be performed with the point of a felt tip pen.
- Execution of Prick-tests**
 - Lay a drop** of each solution to be tested on the indicated spots according to the order on the test sheets.
 - Prick perpendicular to the skin through each extract drop** with the Stallerpoint - one for each test applying a rotation of the device between the thumb and the index. For each prick it is important to ensure that regular and moderate pressure is applied at each time. Moderate pressure is defined as that which produces a depression on the skin between 2 and 5 mm from the bottom of the Stallerpoint. Good execution of the prick is usually observed when a print of the Stallerpoint™ is left on the skin. For best interpretation of the results of the skin reactions, it is advised that each test point be spaced out at about 3 cm.
 - Using the same Stallerpoint to perform several subsequent skin prick tests is associated with false positive results. Use one Stallerpoint for each skin prick test. During a session with a same patient, the user must use the same type of lancet.
- Reading of Prick-tests**
 - Mark on the test sheet the name (Surname and first name) of the patient, add the doctor's seal or stamp, add the date and the time of the prick tests.
 - The reading** is best made, on average, after 20 minutes.
 - Clean the skin** thoroughly with alcohol, taking care not to mix the solutions by sweeping, to remove any mark of the glycerol saline diluents which could impede marking with the felt tip pen.
 - Use a **very fine felt tip pen** to draw an outline of the weal. The presence of the papule can be estimated by touch and the outline can be viewed precisely by use of an oblique light.
 - The test report** is gathered by laying a transparent self-adhesive tape on the skin sketching the outline of the weal down on the test sheets according to the test point.
- 4. Interpretation of results**
 - It is possible when the **skin reactivity** meets two conditions :
 - Positive control response: The weal diameter of a prick test with 9% codeine phosphate (or 10 mg/ml histamine hydrochloride) must be larger than or equal to 3 mm.
 - Negative control response: The weal diameter of a prick test with the glycerol saline diluent must be less than or equal to 3 mm.
 - The test is **positive** when :
 - The weal diameter is higher than 3 mm or 3 mm larger than the negative control one
 - The test is negative when the weal response is less than 3 mm.

Adverse reactions:
Risk of light bleeding at the site. In order to limit this, reduce the applied pressure.
Any serious incident that has occurred in relation to Stallerpoint should be reported to Stallergenes (service.pharmacovigilance@stallergenesreer.com) and the competent authority of the Member State in which the user and/or patient is established.

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