NADAL[®] PROM Amniotic Fluid Test (test strip)

REF 431006N-01/431006N-03/431006N-10/
431006N-20

DE Gebrauchsanweisung	2	CZ Návod k použití	29
EN Instructions for use	6	FI Käyttöohje	32
FR Instructions d'utilisation	10	Symbols	35
ES Instrucciones de uso	14	Our Teams	36
IT Istruzioni per l'uso	17		
PL Sposób użycia	21		
PT Instruções de Utilização	25		

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5. Bell SC, Hales MW, Patel SR, Kirwan PH, Drife JO, Milford-Ward A. Amniotic fluid concentrations of secreted pregnancy-associated endometrial alpha1 and alpha2 globulins (alpha1 and alpha2 PEG). Br J Obstet Gynaecol 1986;93:909-915
6. Suikkari AM, Rutanen EM, Seppälä M. Circulating levels of immunoreactive insulin-like growth factor-binding protein in non-pregnant women. Hum Reprod 1987;2:297-300
7. Bell SC, Keyte JW; N-terminalo amnio acid sequence of human pregnancy-associated endometrial alpha1-globulin, an endometrial insulin-like growth factor (IGF) binding protein – Evidence for two small molecular weight IGF binding proteins. The Endoc Soc 1988:1202-1204
8. Rutanen EM, Pekonen F, Kärkkäinen T. Measurement of insulin-like growth factor binding protein-1 in cervical/vaginal secretions: comparison with the ROM-check membrane immunoassay in the diagnosis of ruptured fetal membranes. Clin Chim Acta 1993;214:73-81
9. Medina TM, Hill DA. Preterm premature rupture of membranes: diagnosis and management. Am Famul Physic 2006;73: 659-664
10. Lockwood CJ, Wein R, Chien D, Ghidini A, Alvarez M, Berkowitz RL. Fetal membrane rupture is associated with the presence of insulin-like growth factor-binding protein-1 in vaginal secretions. Am J Obstet Gynecol 1994;171:146-150
11. Guibourdenche J, Luton D, André E, Noël M, Porquet D. Détection rapide de l'insulin-like growth factor-binding protein-1 et de la fibronectine foetale dans les sécrétions cervico-vaginales pour le diagnostic de la rupture prématurée des membranes. Ann Clin Biochem 1999 ;36 :388-390
12. Sucak A, Moroy P, Cakmak P, Mungan T, Mollamahmutoglu L, Danisman N. Insulin-Like Growth Factor Binding Protein-1; a rapid detection of amniotic fluid leakage after amniocentesis. Turk J Med Sci 2005;35:157-161
13. Talabani H, Recoules A, Grenier S, Luton D, Paper T, Porquet D, Oury JF, Guibourdenche J. L'alpha1-microglobuline placentaire dans les sécrétions vaginales : un nouveau marqueur de la rupture prématurée des membranes. Act Pharm Biol Clin 2005 :12 :233-240
14. Akercan F, Cirpan T, Kazandi M, Terek MC, Mgoyi I, Ozkinay E. The value of insulin-like growth factor binding protein-1 in the cervical-vaginal secretion detected by immunochromatographic dipstick test in the prediction of delivery in women with clinically unconfirmed preterm premature rupture of membranes. Eur J Obstet Gynecol Rep Biol 2005:121:159-163
15. Schuman W. Double sac with secondary rupture of the bag of waters during labor; a clinical entity, and its explanation from examination of the membranes. Am J Obstet Gynecol 1951;62:633-8.
16. Rutanen Eeva-Marja, Kärkkäinen Tytti H., Lehtovirta Johanna, Uotila Jukka T., Hinkula Marianne K., Hartikainen Anna-Liisa. Evaluation of a rapid strip test for insulin-like growth factor binding protein-1 in the diagnosis of ruptured fetal membranes. Clinica Chimica Acta 253 (1996) 91-101.

Rev.0, 2020-11-16 OM/UJa

NADAL® PROM Amniotic Fluid Test

(Ref. 431006N-01/431006N-03/431006N-10/431006N-20)

1. Intended Use

The NADAL® PROM Amniotic Fluid Test is an immunochromatographic rapid test for the qualitative detection of IGFBP-1 (Insulin-like Growth Factor Binding Protein 1), a major protein marker of the amniotic fluid, in a vaginal swab. The test is intended as an aid in the diagnosis of the rupture of fetal membranes in pregnant women and designed for professional use only.

2. Introduction and Clinical Significance

Premature rupture of membranes (PROM) is relatively frequent and occurs in 5 to 10% of all pregnancies. It may lead to premature delivery and fetal infection. Conventional clinical examinations cannot always detect a leakage of amniotic fluid. In some cases, biological tests are required to confirm the suspicion of PROM. Biological tests are based on the detection of alkalization of vaginal secretion (an easy-to-carry-out, sensitive, inexpensive procedure, but with poor specificity) or on the detection of molecules present in high concentrations in amniotic fluid (diamine oxidase, alpha feto protein, fibronectin or IGFBP-1). Detection of PROM with a simple rapid test usable at the point of care enables an immediate decision-making.

Failure to diagnose patients with PROM may result in the failure to implement salutary obstetric measures (e.g. neonatal sepsis occurs if PROM is detected more than 32 hours after initial occurrence). On the other hand, a false diagnosis of PROM can lead to inappropriate interventions (e.g. hospitalization, induction of labour). Therefore, a correct and timely diagnosis of this disorder by the clinician is crucial.


However, accurate diagnosis of PROM remains a challenge for obstetrics. Unfortunately, there is no 'gold standard' for diagnosis of PROM in the clinical practice. Most of the commonly available tests are either inaccurate or to some degree invasive.

The NADAL® PROM Amniotic Fluid Test solves all of these problems. It is a rapid test intended to detect the rupture of fetal membranes, providing highly accurate and timely diagnosis of PROM. Consequently, swift measures can be taken to prevent complications (prophylactic use of antibiotics, tocolytic drugs or corticosteroids, labour induction etc.).

3. Test Principle

The NADAL® PROM Amniotic Fluid Test enables the detection of Insulin-like Growth Factor Binding Protein 1 (IGFBP-1) in vaginal secretion through visual interpretation of colour development on the test strip. Anti-IGFBP-1 antibodies are immobilised in the test line region of the membrane. During the test, the specimen reacts with further anti-IGFBP-1 antibodies conjugated to coloured particles and precoated onto the conjugate pad of the test cassette. The mixture then migrates along the membrane by capillary action and interacts with the reagents on the membrane. If there is sufficient IGFBP-1 in the specimen, a coloured line will form in the test line region of the membrane. The presence of this coloured line indicates a positive result, while its absence indicates a negative result. The appearance of a coloured line in the control line region serves as a procedural control, indicating

that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

Depending on the packaging unit, the NADAL® PROM Amniotic Fluid test kit contains:

- 1/3/10/20 test strips, individually packed with desiccant
- Provided additional material according to 93/42/EEC:
1/3/10/20 sterile Copan flocced swabs (FLOQSwabs™) CE 0123,

 Copan Italia S.p.A.,
Via Perotti 10, 25125 Brescia, Italy

- 1/3/10/20 collection tubes with extraction buffer
- 1 tube holder (except Ref. 431006N-01)
- 1 package insert

5. Additional Materials Required

- Timer

6. Storage & Stability

The NADAL® PROM Amniotic Fluid Test should be stored at 2-30°C until the expiry date printed on the sealed foil pouch. The test strip must remain in the sealed foil pouch until use. Do not freeze tests. Care should be taken to protect components of the test kit from contamination. Do not use the test if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to false results.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.
- Do not use the test beyond the expiration date indicated on the package.
- Do not use the test if the foil pouch is damaged.
- Do not reuse tests.
- The test strip should remain in the sealed foil pouch until use.
- Do not dip the test strip beyond the maximum immersion mark.
- Do not add samples to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area) or the immersion area.

NADAL® PROM Amniotic Fluid Test

(Ref. 431006N-01/431006N-03/431006N-10/431006N-20)

- Avoid cross-contamination of specimens by using a new extraction tube for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not swap caps between different tubes with extraction buffer.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Do not use swabs from damaged pouches.
- Humidity and temperature can adversely affect test results.
- Used testing materials should be discarded according to local regulations.

Leave the test strip in this position.


4. Remove the test strip after 5 minutes. Read the results by placing the test strip on a clean, dry and flat surface. Do not interpret test results after more than 10 minutes.

Note: For a sample collected during a strong leakage of amniotic fluid, a positive result can become visible earlier, while for a sample collected during a very small leakage, it will take the full 10 minutes for a positive result to appear.

10. Result Interpretation**Positive**

Two coloured lines appear on the membrane. One line appears in the control line region (C) and the other line appears in the test line region (T). The result indicates that the concentration of IGFBP-1 is above the detection limit of the test.

**Negative**

A coloured line appears in the control line region (C). No apparent coloured line appears in the test line region (T). The result indicates that the concentration of IGFBP-1 is below the detection limit of the test.

**Invalid**

The control line (C) fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded. Please review the procedure and repeat the test with a new test strip. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

**Note:**

The colour intensity in the test line region (T) may vary depending on the concentration of the analyte present in the specimen. Therefore, any shade of colour in the test line region should be considered positive. Note that this is a qualitative test only and it cannot determine the concentration of the analyte in the specimen.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for the control line failure.

11. Quality Control

An internal procedural control is included in the test strip:

A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Good laboratory practice (GLP) recommends the use of control materials to ensure proper test kit performance.

8. Specimen Collection and Preparation

Use a sterile Copan flocked swab to collect a sample of vaginal secretions from the vaginal wall. Carefully insert a swab into the vagina approximately 5 cm deep for 1 minute. Specimen collection from the cervix should be avoided.

Alternatively, a speculum can be used and vaginal secretion can be collected by placing the swab in contact with the vaginal wall at the level of the posterior fornix for 15 seconds.

9. Test Procedure

Ensure that the test components and samples are at room temperature (15-30°C) prior to testing.

1. Hold the collection tube containing extraction buffer on the lid and shake it thoroughly so that the whole liquid reaches the bottom of the tube. Open the collection tube and place it vertically on a flat, horizontal surface. Use the tube holder as an aid.
2. Dip a swab into the tube and rotate it for 10 seconds. Press the swab against the tube walls in order to squeeze as much liquid as possible from it. Discard the swab.
3. Dip a strip vertically into the tube, ensuring the arrows are pointing downwards, so that it touches the bottom of the tube. This will enhance migration of the sample.

12. Limitations

- The NADAL® PROM Amniotic Fluid Test is for professional *in-vitro* diagnostic use and should be only used for the qualitative detection of IGFBP-1.
- The NADAL® PROM Amniotic Fluid Test only detects the presence of IGFBP-1 in the specimen and should not be used as the sole criterium for the diagnosis of PROM.
- In pregnant women, the IGFBP-1 concentration in the blood is much lower than in the amniotic fluid, meaning that a definitive diagnosis is normally possible. However, if the sample material has been contaminated by blood, false positive results may occur which cannot be attributed to membrane rupture.
- The epithelial cells in the cervix may contain IGFBP-1 concentrations which may cause a false positive test result even without the presence of premature rupture of membranes. This also applies to non-pregnant women, as IGFBP-1 is a major protein in the endometrium around mid-cycle.
- According to existing theories (e.g. 'double sac'), ruptures can occur which may lead to a temporary loss of amniotic fluid. Tests carried out at different times can therefore deliver different results.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended.
- As with all diagnostic tests, a confirmed diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated.

13. Performance Characteristics

Sensitivity

The analytical sensitivity of the NADAL® PROM Amniotic Fluid Test is 5 ng/mL IGFBP-1.

Accuracy

The results obtained using the NADAL® PROM Amniotic Fluid Test were compared to those of another commercially available PROM rapid test in a multi-center clinical evaluation. The results of the study demonstrated 97.7% accuracy of the NADAL® PROM Amniotic Fluid Test when compared to the PROM rapid test of a different manufacturer.

NADAL® PROM Amniotic Fluid Test	PROM rapid test of a different manufacturer		
	Positive	Negative	Total
Positive	83	3	86
Negative	2	133	135
Total	85	136	221

Sensitivity: 97.6%

Specificity: 97.8%

Overall agreement: 97.7%

Interfering substances/microorganisms

The following substances and microorganisms have been evaluated at the concentrations listed below. None of them were found to affect the test performance of the NADAL® PROM Amniotic Fluid Test.

Substance/microorganism	Concentration
Pevaryl®	30 mg/mL
GYNO-TROSYD®	20 mg/mL
Flagyl®	100 mg/mL
Canesten®	40 mg/dL
Glucose	8 mg/mL
<i>Candida albicans</i>	11.2 x 10 ⁸ CFU/mL
<i>Gardnerella vaginalis</i>	8.6 x 10 ⁸ CFU/mL
<i>Neisseria gonorrhoeae</i>	10.6 x 10 ⁸ CFU/mL
Baby powder	10 mg/mL
Vagisan® Myco Kombi (Creme)	1 mg/mL
Baby oil	1 mg/mL
Medical silicone oil	13 g/mL

Note: blood of pregnant women may cause false positive results.

Intra-lot and inter-lot variability

Negative and positive controls at concentrations of 5 ng/mL and 25 ng/mL were tested in 10 replicates with 3 independent lots on 3 consecutive days. Each control was recognized correctly.

14. References

1. Friedman ML, McElin TW. Diagnosis of ruptured fetal membranes. Am J Obstet Gynecol 1969;104:544-550
2. Steinman G, Kleiner GJ, Greston WM. Spontaneous rupture of membranes. NY St. J Med 1979;1849-1851
3. Koninckx PR, Trappeniens H, van Assche FA. Prolactin concentration in vaginal fluid: a new method for diagnosing ruptured membranes. Br J Obstet Gynaecol 1981;88:607-610
4. Rochelson BL, Richardson DA, Macri JN. Rapid assay: possible application in the diagnosis of premature rupture of the membranes. Obstet Gynecol 1983;62:414-418
5. Bell SC, Hales MW, Patel SR, Kirwan PH, Drife JO, Milford-Ward A. Amniotic fluid concentrations of secreted pregnancy-associated endometrial alpha1 and alpha2 globulins (alpha1 and alpha2 PEG). Br J Obstet Gynaecol 1986;93:909-915
6. Suikari AM, Rutanen EM, Seppälä M. Circulating levels of immunoreactive insulin-like growth factor-binding protein in non-pregnant women. Hum Reprod 1987;2:297-300
7. Bell SC, Keyte JW; N-terminal amino acid sequence of human pregnancy-associated endometrial alpha1-globulin, an endometrial insulin-like growth factor (IGF) binding protein - Evidence for two small molecular weight IGF binding proteins. The Endoc Soc 1988:1202-1204
8. Rutanen EM, Pekonen F, Kärkkäinen T. Measurement of insulin-like growth factor binding protein-1 in cervical/vaginal secretions: comparison with the ROM-check membrane immunoassay in the diagnosis of ruptured fetal membranes. Clin Chim Acta 1993;214:73-81
9. Medina TM, Hill DA. Preterm premature rupture of membranes: diagnosis and management. Am Fam Physic 2006;73: 659-664
10. Lockwood CJ, Wein R, Chien D, Ghidini A, Alvarez M, Berkowitz RL. Fetal membrane rupture is associated with the presence of insulin-like growth factor-binding protein-1 in vaginal secretions. Am J Obstet Gynecol 1994;171:146-150
11. Guibourdenche J, Luton D, André E, Noël M, Porquet D. Détection rapide de l'insulin-like growth factor-binding protein-1 et de la fibronectine fœtale dans les sécrétions cervico-vaginales pour le diagnostic de la rupture prématurée des membranes. Ann Clin Biochem 1999;36:388-390
12. Suckak A, Moroy P, Cakmak P, Mungan T, Mollamahmutoglu L, Danisman N. Insulin-Like Growth Factor Binding Protein-I: a rapid detection of amniotic fluid leakage after amniocentesis. Turk J Med Sci 2005;35:157-161
13. Talabani H, Recoules A, Grenier S, Luton D, Paper T, Porquet D, Oury JF, Guibourdenche J. L'alpha1-microglobuline placentaire dans les sécrétions vaginales : un nouveau marqueur de la rupture prématurée des membranes. Act Pharm Biol Clin 2005 ;12 :233-240
14. Akercan F, Cirpan T, Kazandi M, Terek MC, Migoyi I, Ozkinay E. The value of insulin-like growth factor binding protein-1 in the cervical-vaginal secretion detected by immunochromatographic dipstick test in the prediction of delivery in women with clinically unconfirmed preterm premature rupture of membranes. Eur J Obstet Gynecol Rep Biol 2005;121:159-163
15. Schuman W. Double sac with secondary rupture of the bag of waters during labor; a clinical entity, and its explanation from examination of the membranes. Am J Obstet Gynecol 1951;62:633-8.

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7. Bell SC, Keyte JW; N-terminalo amnio acid sequence of human pregnancy-associated endometrial alpha1-globulin, an endometrial insulin-like growth factor (IGF) binding protein – Evidence for two small molecular weight IGF binding proteins. *The Endoc Soc* 1988;1202-1204
8. Rutanen EM, Pekonen F, Kärkkäinen T. Measurement of insulin-like growth factor binding protein-1 in cervical/vaginal secretions: comparison with the ROM-check membrane immunoassay in the diagnosis of ruptured fetal membranes. *Clin Chim Acta* 1993;214:73-81
9. Medina TM, Hill DA. Preterm premature rupture of membranes: diagnosis and management. *Am Fam Physic* 2006;73: 659-664
10. Lockwood CJ, Wein R, Chien D, Ghidini A, Alvarez M, Berkowitz RL. Fetal membrane rupture is associated with the presence of insulin-like growth factor-binding protein-1 in vaginal secretions. *Am J Obstet Gynecol* 1994;171:146–150
11. Guibourdenche J, Luton D, André E, Noël M, Porquet D. Détection rapide de l'insulin-like growth factor-binding protein-1 et de la fibronectine fœtale dans les sécrétions cervico-vaginales pour le diagnostic de la rupture prématurée des membranes. *Ann Clin Biochem* 1999 ;36 :388-390
12. Sucak A, Moroy P, Cakmak P, Mungan T, Mollamahmutoglu L, Danisman N. Insulin-Like Growth Factor Binding Protein-1; a rapid detection of amniotic fluid leakage after amniocentesis. *Turk J Med Sci* 2005;35:157-161
13. Talabani H, Recoules A, Grenier S, Luton D, Paper T, Porquet D, Oury JF, Guibourdenche J. L'alpha1-microglobuline placentaire dans les sécrétions vaginales : un nouveau marqueur de la rupture prématurée des membranes. *Act Pharm Biol Clin* 2005 ;12 :233-240
14. Akercan F, Cirpan T, Kazandi M, Terek MC, Mgoyl I, Ozkinay E. The value of insulin-like growth factor binding protein-1 in the cervical-vaginal secretion detected by immunochromatographic dipstick test in the prediction of delivery in women with clinically unconfirmed preterm premature rupture of membranes. *Eur J Obstet Gynecol Rep Biol* 2005;121:159-163
15. Schuman W. Double sac with secondary rupture of the bag of waters during labor; a clinical entity, and its explanation from examination of the membranes. *Am J Obstet Gynecol* 1951;62:633–8.
16. Rutanen Eeva-Marja, Kärkkäinen Tytti H., Lehtovirta Johanna, Uotila Jukka T., Hinkula Marianne K., Hartikainen Anna-Liisa. Evaluation of a rapid strip test for insulin-like growth factor binding protein-1 in the diagnosis of ruptured fetal membranes. *Clinica Chimica Acta* 253 (1996) 91-101.

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14. Referencias

- Friedman ML, McElin TW. Diagnosis of ruptured fetal membranes. Am J Obstet Gynecol 1969;104:544-550
- Steinman G, Kleiner GJ, Greston WM. Spontaneous rupture of membranes. NY St. J Med 1979;1849-1851
- Koninckx PR, Trappeniers H, van Assche FA. Prolactin concentration in vaginal fluid: a new method for diagnosing ruptured membranes. Br J Obstet Gynaecol 1981;88:607-610
- Rochelson BL, Richardson DA, Macri JN. Rapid assay: possible application in the diagnosis of premature rupture of the membranes. Obstet Gynecol 1983;62:414-418
- Bell SC, Hales MW, Patel SR, Kirwan PH, Drife JO, Milford-Ward A. Amniotic fluid concentrations of secreted pregnancy-associated endometrial alpha1 and alpha2 globulins (alpha1 and alpha2 PEG). Br J Obstet Gynaecol 1986;93:909-915
- Suikkari AM, Rutanen EM, Seppälä M. Circulating levels of immunoreactive insulin-like growth factor-binding protein in non-pregnant women. Hum Reprod 1987;2:297-300
- Bell SC, Keyte JW; N-terminal amino acid sequence of human pregnancy-associated endometrial alpha1-globulin, an endometrial insulin-like growth factor (IGF) binding protein - Evidence for two small molecular weight IGF binding proteins. The Endoc Soc 1988:1202-1204
- Rutanen EM, Pekonen F, Kärkkäinen T. Measurement of insulin-like growth factor binding protein-1 in cervical/vaginal secretions: comparison with the ROM-check membrane immunoassay in the diagnosis of ruptured fetal membranes. Clin Chim Acta 1993;214:73-81
- Medina TM, Hill DA. Preterm premature rupture of membranes: diagnosis and management. Am Fam Physic 2006;73: 659-664
- Lockwood CJ, Wein R, Chien D, Ghidini A, Alvarez M, Berkowitz RL. Fetal membrane rupture is associated with the presence of insulin-like growth factor-binding protein-1 in vaginal secretions. Am J Obstet Gynecol 1994;171:146-150
- Guibourdenche J, Luton D, André E, Noël M, Porquet D. Détection rapide de l'insuline-like growth factor-binding protein-1 et de la fibronectine fœtale dans les sécrétions cervico-vaginales pour le diagnostic de la rupture prématurée des membranes. Ann Clin Biochem 1999;36:388-390
- Sucak A, Moroy P, Cakmak P, Mungan T, Mollamahmutoglu I, Danisman N. Insulin-Like Growth Factor Binding Protein-1: a rapid detection of amniotic fluid leakage after amniocentesis. Turk J Med Sci 2005;35:157-161
- Talabani H, Recoules A, Grenier S, Luton D, Paper T, Porquet D, Oury JF, Guibourdenche J. L'alpha1-microglobuline placentaire dans les sécrétions vaginales: un nouveau marqueur de la rupture prématurée des membranes. Act Pharm Biol Clin 2005;12:233-240
- Akercan F, Cirpan T, Kazandi M, Terek MC, Migoyi I, Ozkinay E. The value of insulin-like growth factor binding protein-1 in the cervical-vaginal secretion detected by immunochromatographic dipstick test in the prediction of delivery in women with clinically unconfirmed preterm premature rupture of membranes. Eur J Obstet Gynecol Rep Biol 2005;121:159-163
- Schuman W. Double sac with secondary rupture of the bag of waters during labor; a clinical entity, and its explanation from examination of the membranes. Am J Obstet Gynecol 1951;62:633-8.
- Rutanen Eeva-Marja, Kärkkäinen Tytti H., Lehtovirta Johanna, Uotila Jukka T., Hinkula Marianne K., Hartikainen Anna-Liisa. Evaluation of a rapid strip test for insulin-like growth factor binding protein-1 in the diagnosis of ruptured fetal membranes. Clinica Chimica Acta 253 (1996) 91-101.

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14. Bibliografia

- Friedman ML, McElin TW. Diagnosis of ruptured fetal membranes. Am J Obstet Gynecol 1969;104:544-550
- Steinman G, Kleiner GJ, Greston WM. Spontaneous rupture of membranes. NY St. J Med 1979;1849-1851
- Koninckx PR, Trappeniers H, van Assche FA. Prolactin concentration in vaginal fluid: a new method for diagnosing ruptured membranes. Br J Obstet Gynaecol 1981;88:607-610
- Rochelson BL, Richardson DA, Macri JN. Rapid assay: possible application in the diagnosis of premature rupture of the membranes. Obstet Gynecol 1983;62:414-418
- Bell SC, Hales MW, Patel SR, Kirwan PH, Drife JO, Milford-Ward A. Amniotic fluid concentrations of secreted pregnancy-associated endometrial alpha1 and alpha2 globulins (alpha1 and alpha2 PEG). Br J Obstet Gynaecol 1986;93:909-915
- Suikkari AM, Rutanen EM, Seppälä M. Circulating levels of immunoreactive insulin-like growth factor-binding protein in non-pregnant women. Hum Reprod 1987;2:297-300
- Bell SC, Keyte JW. N-terminal amino acid sequence of human pregnancy-associated endometrial alpha1-globulin, an endometrial insulin-like growth factor (IGF) binding protein - Evidence for two small molecular weight IGF binding proteins. The Endoc Soc 1988;1202-1204
- Rutanen EM, Pekonen F, Kärkkäinen T. Measurement of insulin-like growth factor binding protein-1 in cervical/vaginal secretions: comparison with the ROM-check membrane immunoassay in the diagnosis of ruptured fetal membranes. Clin Chim Acta 1993;214:73-81
- Medina TM, Hill DA. Preterm premature rupture of membranes: diagnosis and management. Am Famul Physic 2006;73: 659-664
- Lockwood CJ, Wein R, Chien D, Ghidini A, Alvarez M, Berkowitz RL. Fetal membrane rupture is associated with the presence of insulin-like growth factor-binding protein-1 in vaginal secretions. Am J Obstet Gynecol 1994;171:146-150
- Guibourdenche J, Luton D, André E, Noël M, Porquet D. Détection rapide de l'Insulin-like growth factor-binding protein-1 et de la fibronectine fœtale dans les sécrétions cervico-vaginales pour le diagnostic de la rupture prématurée des membranes. Ann Clin Biochem 1999;36:388-390
- Sukac A, Moroy P, Cakmak P, Mungan T, Mollamahmutoglu L, Danisman N. Insulin-Like Growth Factor Binding Protein-I; a rapid detection of amniotic fluid leakage after amniocentesis. Turk J Med Sci 2005;35:157-161

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14. Bibliografia

1. Friedman ML, McEllin TW. Diagnosis of ruptured fetal membranes. Am J Obstet Gynecol 1969;104:544-550
2. Steinman G, Kleiner GJ, Greston WM. Spontaneous rupture of membranes. NY St. J Med 1979;1849-1851
3. Koninckx PR, Trappeniers H, van Assche FA. Prolactin concentration in vaginafluid: a new method for diagnosing ruptured membranes. Br J Obstet Gynaecol 1981;88:607-610
4. Rochelson BL, Richardson DA, Macri JN. Rapid assay: possible application in the diagnosis of premature rupture of the membranes. Obstet Gynecol 1983;62:414-418
5. Bell SC, Hales MW, Patel SR, Kirwan PH, Drife JO, Milford-Ward A. Amniotic fluid concentrations of secreted pregnancy-associated endometrial alpha1 and alpha2 globulins (alpha1 and alpha2 PEG). Br J Obstet Gynaecol 1986;93:909-915
6. Suikari AM, Rutanen EM, Seppälä M. Circulating levels of immunoreactive insulin-like growth factor-binding protein in non-pregnant women. Hum Reprod 1987;2:297-300
7. Bell SC, Kayte JW; N-terminal amino acid sequence of human pregnancy-associated endometrial alpha1-globulin, an endometrial insulin-like growth factor (IGF) binding protein – Evidence for two small molecular weight IGF binding proteins. The Endoc Soc 1988:1202-1204
8. Rutanen EM, Pekonen F, Kärkkäinen T. Measurement of insulin-like growth factor binding protein-1 in cervical/vaginal secretions: comparison with the ROM-check

- membrane immunoassay in the diagnosis of ruptured fetal membranes. Clin Chim Acta 1993;214:73-81
9. Medina TM, Hill DA. Preterm premature rupture of membranes: diagnosis and management. Am Fam Physic 2006;73: 659-664
 10. Lockwood CJ, Wein R, Chien D, Ghidini A, Alvarez M, Berkowitz RL. Fetal membrane rupture is associated with the presence of insulin-like growth factor-binding protein-1 in vaginal secretions. Am J Obstet Gynecol 1994;171:146-150
 11. Guibourdenche J, Luton D, André E, Noël M, Porquet D. Détection rapide de l'insulin-like growth factor-binding protein-1 et de la fibronectine fœtale dans les sécrétions cervico-vaginales pour le diagnostic de la rupture prématurée des membranes. Ann Clin Biochem 1999 ;36 :388-390
 12. Sucak A, Moroy P, Cakmak P, Mungan T, Mollamahmutoglu L, Danisman N. Insulin-Like Growth Factor Binding Protein-1; a rapid detection of amniotic fluid leakage after amniocentesis. Turk J Med Sci 2005;35:157-161
 13. Talabani H, Recoules A, Grenier S, Luton D, Paper T, Porquet D, Oury JF, Guibourdenche J. L'alpha1-microglobuline placentaire dans les sécrétions vaginales : un nouveau marqueur de la rupture prématurée des membranes. Act Pharm Biol Clin 2005 :12 :233-240
 14. Akercan F, Cirpan T, Kazandi M, Terek MC, Mgoyi I, Ozkinay E. The value of insulin-like growth factor binding protein-1 in the cervical-vaginal secretion detected by immunochromatographic dipstick test in the prediction of delivery in women with clinically unconfirmed preterm premature rupture of membranes. Eur J Obstet Gynecol Rep Biol 2005;121:159-163
 15. Schuman W. Double sac with secondary rupture of the bag of waters during labor; a clinical entity, and its explanation from examination of the membranes. Am J Obstet Gynecol 1951;62:633-8.
 16. Rutanen Eeva-Marja, Kärkkäinen Tytti H., Lehtovirta Johanna, Uotila Jukka T., Hinkula Marianne K., Hartikainen Anna-Liisa. Evaluation of a rapid strip test for insulin-like growth factor binding protein-1 in the diagnosis of ruptured fetal membranes. Clinica Chimica Acta 253 (1996) 91-101.

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14. Referências

1. Friedman ML, McEllin TW. Diagnosis of ruptured fetal membranes. Am J Obstet Gynecol 1969;104:544-550
2. Steinman G, Kleiner GJ, Greston WM. Spontaneous rupture of membranes. NY St. J Med 1979;1849-1851
3. Koninckx PR, Trappeniers H, van Assche FA. Prolactin concentration in vaginafluid: a new method for diagnosing ruptured membranes. Br J Obstet Gynaecol 1981;88:607-610
4. Rochelson BL, Richardson DA, Macri JN. Rapid assay: possible application in the diagnosis of premature rupture of the membranes. Obstet Gynecol 1983;62:414-418
5. Bell SC, Hales MW, Patel SR, Kirwan PH, Drife JO, Milford-Ward A. Amniotic fluid concentrations of secreted pregnancy-associated endometrial alpha1 and alpha2 globulins (alpha1 and alpha2 PEG). Br J Obstet Gynaecol 1986;93:909-915
6. Suikkari AM, Rutanen EM, Seppälä M. Circulating levels of immunoreactive insulin-like growth factor-binding protein in non-pregnant women. Hum Reprod 1987;2:297-300
7. Bell SC, Keyte JW; N-terminal amino acid sequence of human pregnancy-associated endometrial alpha1-globulin, an endometrial insulin-like growth factor (IGF) binding protein - Evidence for two small molecular weight IGF binding proteins. The Endoc Soc 1988:1202-1204
8. Rutanen EM, Pekonen F, Kärkäinen T. Measurement of insulin-like growth factor binding protein-1 in cervical/vaginal secretions: comparison with the ROM-check membrane immunoassay in the diagnosis of ruptured fetal membranes. Clin Chim Acta 1993;214:73-81
9. Medina TM, Hill DA. Preterm premature rupture of membranes: diagnosis and management. Am Fam Physic 2006;73: 659-664

10. Lockwood CJ, Wein R, Chien D, Ghidini A, Alvarez M, Berkowitz RL. Fetal membrane rupture is associated with the presence of insulin-like growth factor-binding protein-1 in vaginal secretions. *Am J Obstet Gynecol* 1994;171:146-150
11. Guibourdenche J, Luton D, André E, Noël M, Porquet D. Détection rapide de l'insulin-like growth factor-binding protein-1 et de la fibronectine fœtale dans les sécrétions cervico-vaginales pour le diagnostic de la rupture prématurée des membranes. *Ann Clin Biochem* 1999 ;36 :388-390
12. Sucak A, Moroy P, Cakmak P, Mungan T, Mollamahmutoglu L, Danisman N. Insulin-Like Growth Factor Binding Protein-1; a rapid detection of amniotic fluid leakage after amniocentesis. *Turk J Med Sci* 2005;35:157-161
13. Talabani H, Recoules A, Grenier S, Luton D, Paper T, Porquet D, Oury JF, Guibourdenche J. L'alpha1-microglobuline placentaire dans les sécrétions vaginales : un nouveau marqueur de la rupture prématurée des membranes. *Act Pharm Biol Clin* 2005 :12 :233-240
14. Akercan F, Cirpan T, Kazandi M, Terek MC, Mgoyi I, Ozkinay E. The value of insulin-like growth factor binding protein-1 in the cervical-vaginal secretion detected by immunochromatographic dipstick test in the prediction of delivery in women with clinically unconfirmed preterm premature rupture of membranes. *Eur J Obstet Gynecol Rep Biol* 2005;121:159-163
15. Schuman W. Double sac with secondary rupture of the bag of waters during labor; a clinical entity, and its explanation from examination of the membranes. *Am J Obstet Gynecol* 1951;62:633-8.
16. Rutanen Eeva-Marja, Kärkkäinen Tytti H., Lehtovirta Johanna, Uotila Jukka T., Hinkula Marianne K., Hartikainen Anna-Liisa. Evaluation of a rapid strip test for insulin-like growth factor binding protein-1 in the diagnosis of ruptured fetal membranes. *Clinica Chimica Acta* 253 (1996) 91-101.

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14. Reference

- Friedman ML, McElin TW. Diagnosis of ruptured fetal membranes. Am J Obstet Gynecol 1969;104:544-550
- Steinman G, Kleiner GJ, Greston WM. Spontaneous rupture of membranes. NY St. J Med 1979;1849-1851
- Koninckx PR, Trappeniers H, van Assche FA. Prolactin concentration in vaginal fluid: a new method for diagnosing ruptured membranes. Br J Obstet Gynaecol 1981;88:607-610
- Rochelson BL, Richardson DA, Macri JN. Rapid assay: possible application in the diagnosis of premature rupture of the membranes. Obstet Gynecol 1983;62:414-418
- Bell SC, Hales MW, Patel SR, Kirwan PH, Drife JO, Milford-Ward A. Amniotic fluid concentrations of secreted pregnancy-associated endometrial alpha1 and alpha2 globulins (alpha1 and alpha2 PEG). Br J Obstet Gynaecol 1986;93:909-915
- Suikkari AM, Rutanen EM, Seppälä M. Circulating levels of immunoreactive insulin-like growth factor-binding protein in non-pregnant women. Hum Reprod 1987;2:297-300
- Bell SC, Keyte JW; N-terminal amnio acid sequence of human pregnancy-associated endometrial alpha1-globulin, an endometrial insulin-like growth factor (IGF) binding protein - Evidence for two small molecular weight IGF binding proteins. The Endoc Soc 1988:1202-1204
- Rutanen EM, Pekonen E, Kärkkäinen T. Measurement of insulin-like growth factor binding protein-1 in cervical/vaginal secretions: comparison with the ROM-check membrane immunoassay in the diagnosis of ruptured fetal membranes. Clin Chim Acta 1993;214:73-81
- Medina TM, Hill DA. Preterm premature rupture of membranes: diagnosis and management. Am Famul Physic 2006;73: 659-664
- Lockwood CJ, Wein R, Chien D, Ghidini A, Alvarez M, Berkowitz RL. Fetal membrane rupture is associated with the presence of insulin-like growth factor-binding protein-1 in vaginal secretions. Am J Obstet Gynecol 1994;171:146-150
- Guibourdenne J, Luton D, André E, Noël M, Porquet D. Détection rapide de l'insulin-like growth factor-binding protein-1 et de la fibronectine fœtale dans les sécrétions cervico-vaginales pour le diagnostic de la rupture prématurée des membranes. Ann Clin Biochem 1999 ;36 :388-390
- Sucak A, Moroy P, Cakmak P, Mungan T, Mollamahmutoglu L, Danisman N. Insulin-Like Growth Factor Binding Protein-I; a rapid detection of amniotic fluid leakage after amniocentesis. Turk J Med Sci 2005;35:157-161
- Talabani H, Recoules A, Grenier S, Luton D, Papper T, Porquet D, Oury JF, Guibourdenne J. L'alpha1-microglobuline placentaire dans les sécrétions vaginales : un nouveau marqueur de la rupture prématurée des membranes. Act Pharm Biol Clin 2005 ;12 :233-240
- Akercan F, Cirpan T, Kazandi M, Terek MC, Mgoiy I, Ozkinay E. The value of insulin-like growth factor binding protein-1 in the cervical-vaginal secretion detected by immunochromatographic dipstick test in the prediction of delivery in women with clinically unconfirmed preterm premature rupture of membranes. Eur J Obstet Gynecol Rep Biol 2005;121:159-163
- Schuman W. Double sac with secondary rupture of the bag of waters during labor; a clinical entity, and its explanation from examination of the membranes. Am J Obstet Gynecol 1951;62:633-8.
- Rutanen Eeva-Marja, Kärkkäinen Tytti H., Lehtovirta Johanna, Uotila Jukka T., Hinkula Marianne K., Hartikainen Anna-Liisa. Evaluation of a rapid strip test for insulin-like growth factor binding protein-1 in the diagnosis of ruptured fetal membranes. Clinica Chimica Acta 253 (1996) 91-101.

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



14. Lähteet

1. Friedman ML, McElin TW. Diagnosis of ruptured fetal membranes. Am J Obstet Gynecol 1969;104:544-550
2. Steinman G, Kleiner GJ, Grestone WM. Spontaneous rupture of membranes. NY St. J Med 1979;1849-1851
3. Koninckx PR, Trappeniers H, van Assche FA. Prolactin concentration in vaginal fluid: a new method for diagnosing ruptured membranes. Br J Obstet Gynaecol 1981;88:607-610
4. Rochelson BL, Richardson DA, Macri JN. Rapid assay: possible application in the diagnosis of premature rupture of the membranes. Obstet Gynecol 1983;62:414-418
5. Bell SC, Hales MW, Patel SR, Kirwan PH, Drife JO, Milford-Ward A. Amniotic fluid concentrations of secreted pregnancy-associated endometrial alpha1 and alpha2 globulins (alpha1 and alpha2 PEG). Br J Obstet Gynaecol 1986;93:909-915
6. Suikkari AM, Rutanen EM, Seppälä M. Circulating levels of immunoreactive insulin-like growth factor-binding protein in non-pregnant women. Hum Reprod 1987;2:297-300
7. Bell SC, Keyte JW; N-terminalo amnio acid sequence of human pregnancy-associated endometrial alpha1-globulin, an endometrial insulin-like growth factor (IGF) binding protein - Evidence for two small molecular weight IGF binding proteins. The Endoc Soc 1988:1202-1204
8. Rutanen EM, Pekonen F, Kärkkäinen T. Measurement of insulin-like growth factor binding protein-1 in cervical/vaginal secretions: comparison with the ROM-check membrane immunoassay in the diagnosis of ruptured fetal membranes. Clin Chim Acta 1993;214:73-81
9. Medina TM, Hill DA. Preterm premature rupture of membranes: diagnosis and management. Am Fam Physic 2006;73: 659-664
10. Lockwood CJ, Wein R, Chien D, Ghidini A, Alvarez M, Berkowitz RL. Fetal membrane rupture is associated with the presence of insulin-like growth factor-binding protein-1 in vaginal secretions. Am J Obstet Gynecol 1994;171:146-150
11. Guibourdenche J, Luton D, André E, Noël M, Porquet D. Détection rapide de l'Insulin-like growth factor-binding protein-1 et de la fibronectine fœtale dans les sécrétions cervico-vaginales pour le diagnostic de la rupture prématurée des membranes. Ann Clin Biochem 1999 ;36 :388-390
12. Suckak A, Moroy P, Cakmak P, Mungan T, Mollamahmutoglu L, Danisman N. Insulin-Like Growth Factor Binding Protein-I: a rapid detection of amniotic fluid leakage after amniocentesis. Turk J Med Sci 2005;35:157-161
13. Talabani H, Recoules A, Grenier S, Luton D, Papper T, Porquet D, Oury JF, Guibourdenche J. L'alpha1-microglobuline placentaire dans les sécrétions vaginales : un nouveau marqueur de la rupture prématurée des membranes. Act Pharm Biol Clin 2005 ;12 :233-240
14. Akercan F, Girpan T, Kazandi M, Terek MC, Mgoyi I, Ozkinay E. The value of insulin-like growth factor binding protein-1 in the cervical-vaginal secretion detected by immunochromatographic dipstick test in the prediction of delivery in women with clinically unconfirmed preterm premature rupture of membranes. Eur J Obstet Gynecol Rep Biol 2005;121:159-163
15. Schuman W. Double sac with secondary rupture of the bag of waters during labor: a clinical entity, and its explanation from examination of the membranes. Am J Obstet Gynecol 1951;62:633-8.
16. Rutanen Eeva-Marja, Kärkkäinen Tytti H., Lehtovirta Johanna, Uotila Jukka T., Hinkula Marianne K., Hartikainen Anna-Liisa. Evaluation of a rapid strip test for insulin-like growth factor binding protein-1 in the diagnosis of ruptured fetal membranes. Clinica Chimica Acta 253 (1996) 91-101.

Rev. 0, 2020-11-16 OL

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conforme aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consultense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	<i>in-vitro</i> -Diagnostika	<i>in-vitro</i> diagnostic medical device	Dispositif médical de diagnostic <i>in-vitro</i>	Producto sanitario para diagnóstico <i>in-vitro</i>	Dispositivo medico-diagnostico <i>in-vitro</i>	Tylko do diagnostyki <i>in-vitro</i>
	Temperaturbegrenzung	Temperature limitation	Limites de température	Límite de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Numéro de lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour "n" tests	Suficiente para <n> utilizaciones	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń

Symbol	Português	Český	Suomi	Svenskt	Nederlands	Dansk	Norsk
	Conformidade com as normas europeias	CE certifikát	CE-merkitty	CE-märkning	CE-markering	CE-mærkning	CE standardisert
	Consultar as instruções de utilização	Viz návod k použití	Katso käyttöohjetta	Läs bruksanvisningen	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen	Les bruksanvisning nøye
	Dispositivo médico para diagnóstico <i>in-vitro</i>	Diagnostický zdravotnický prostředek <i>in-vitro</i>	<i>in-vitro</i> - diagnostiikkaan tarkoitettu lääkinnällinen laite	Medicinteknisk produkt avsedd för <i>in-vitro</i> -diagnostik	Medisch hulpmiddel voor <i>in-vitro</i> -diagnostiek	Medicinsk udstyr til <i>in-vitro</i> -diagnostik	<i>in-vitro</i> diagnostic medisinsk enhet
	Limites de temperatura	Teplotní omezení	Lämpötilarajat	Temperaturbegränsning	Temperatuurlimiet	Temperaturbegrænsning	Temperaturbegrensning
	Código do lote	Kód šarže	Eräkoodi	Satsnummer	Code van de partij	Batchkode	Merkning
	Não reutilizar	Pro jednorázové použití	Kertakäyttöinen	Får inte återanvändas	Niet opnieuw gebruiken	Må ikke genbruges	Må ikke brukes om igjen
	Prazo de validade	Spotřebujte do	Käytettävä viimeistään	Används före	Houdbaar tot	Udløbsdato	Tidtaking
	Número de catálogo	Katalogov číslo	Luettelonumero	Listnummer	Catalogus nummer	Best il lingsnummer	Katalog nummer
	Fabricante	Výrobce	Valmistaja	Tillverkare	Fabrikant	Fabrikant	Producent
	Suficiente para <n> test	Dostačuje pro <n> testů	Lukumäärä <n> test	Räcker till <n> test	Voldoende voor <n> test	Tilstrækkeligt til <n> test	Tilstrækkelig for<n> tester

Our Teams

Germany:

Regensburg

Tel: +49 941 290 10-0

Fax: +49 941 290 10-50

Moers

Tel: +49 2841 99820-0

Fax: +49 2841 99820-1

Austria:

Tel: +49 941 290 10-29

Free Tel: 0800 291 565

Fax: +49 290 10-50

Free Fax: 0800 298 197

UK & Ireland:

Tel: +49 941 290 10-18

Free Tel – UK: 0808 234 1237

Free Tel – IRE: 1800 555 080

Fax: +49 290 10-50

France:

France Tel: 0800 915 240

France Fax: 0800 909 493

Switzerland

Swiss Tel: 0800 564 720

Swiss Fax: 0800 837 476

Belgium

Belgium Tel: 0800 718 82

Belgium Fax: 0800 747 07

Luxembourg

Lux, Tel: 800 211 16

Lux, Fax: 800 261 79

Spain:

Tel: +49 941 290 10-759

Free Tel: 900 938 315

Fax: +49 941 290 10-50

Free Fax: 900 984 992

Italy:

Tel: +49 941 290 10-34

Fax: +49 941 290 10-50

Poland:

Tel: +49 941 290 10-44

Free Tel: 00 800 491 15 95

Fax: +49 941 290 10-50

Free Fax: 00 800 491 15 94

Portugal:

Tel: +49 941 290 10-735

Tel, Verde: 800 849 230

Fax: +49 941 290 10-50

Fax Verde: 800 849 229

Netherlands:

Tel: +31 30 75 600

Free Tel: 0800 0222 890

Fax: +31 70 30 30 775

Free Fax: 0800 024 9519

Nordic countries:

Denmark

Tel: +31 703075 605

Free Tel: 808 887 53

Finland

Tel: +31 703075 606

Free Tel: 0800 918 263

Free Fax: 0800 918 262

Norway

Tel: +31 703075 605

Free Tel: 800 16 731

Sweden

Tel: +31 703075 605

Free Tel: 020 79 09 06

Laboratory Diagnostics Team:

Tel: +49 941 290 10-40

Fax: +49 941 290 10-50

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany

www.nal-vonminden.com • info@nal-vonminden.com

Fon: +49 2841 99820-0 • Fax: +49 2841 99820-1

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nal von minden GmbH