
Introcan Safety[®] 3

B | BRAUN


LLDorder 5590 - grunewald 211068

GB	Closed IV Catheter	14
DE	Geschlossener IV-Katheter	16
BG	Затворен интравенозен катетър	19
BR	Cateter intravascular fechado	22
CN	闭合 IV 导管	24
CZ	Uzavřený intravenózní katétr	26
DK	Lukket IV-kateter	29
EE	Suletud IV-kateeter	31
ES	Catéter intravascular cerrado	34
FI	Suljettu IV-katetri	37
FR	Cathéter IV fermé	39
GR	ΕΦ καθετήρας κλειστού τύπου	42
HR	Zatvoreni intravenski kateter	45
HU	Zárt intravénás katéter	48
IT	Catetere intravascolare a sistema chiuso	51
KR	폐쇄형 IV 카테터	53
LT	Uždaras IV kateteris	56
LV	Slēgts intravenozais katetrs	58
NL	Gesloten IV-katheter	61
NO	Lukket IV-kateter	64
PL	Zamknięty cewnik IV	66

PT	Cateter Intravenoso Fechado	69
RO	Cateter IV închis	71
RS	Zatvoreni IV kateter	74
RU	Безопасный в/в катетер	76
SE	Sluten IV-kateter	79
SI	Zaprta kateter IV	82
SK	Uzavretý katéter IV	84
TH	สายสวนหลอดเลือดดำระบบปิด	87
TR	Kapalı IV Kateter	89
VN	Ống thông Nội mạch Kín	92

GB Fig. 1: Before needle removal
DE Abb. 1: Vor dem Entfernen der Kanüle
BG Фиг. 1: Преди изваждане на иглата
BR Fig. 1: Antes da remoção da agulha
CN 图 1: 撤出针芯之前
CZ Obr. 1: Před vytažením jehly
DK Fig. 1: Inden fjernelse af kanylen
EE Joon. 1: Enne nõela eemaldamist
ES Fig. 1: antes de la extracción de la aguja
FI Kuva 1: Ennen neulan poistamista
FR Fig. 1 : Avant le retrait de l'aiguille
GR Σχ. 1: Πριν από την αφαίρεση της βελόνας
HR Slika 1: prije vađenja igle
HU 1. ábra: A tű eltávolítása előtt
IT Fig. 1: Prima della rimozione dell'ago
KR 그림 1: 바늘 제거 전
LT 1 pav.: prieš nuimant adatą
LV 1. att. Pirms adatas noņemšanas
NL Afb. 1: Voor verwijderen van de naald
NO Fig. 1: Før fjerning av nål
PL Ryc. 1: Przed wyjęciem igły
PT Fig. 1: Antes da remoção da agulha
RO Fig. 1: înainte de scoaterea acului
RS Sl. 1: Pre uklanjanja igle

RU Рис. 1. До извлечения иглы
SE Fig. 1: Innan nålen tas bort
SI Slika 1: Pred odstranitvijo igle
SK Obr. 1: Pred vytiahnutím ihly
TH รูปที่ 1 ก่อนถอดเข็ม
TR Şek. 1: İğneyi çıkarmadan önce
VN Hình 1: Trước khi tháo kim

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





GB Instructions for use

Description:

The Introcath Safety® 3 Closed IV Catheter is designed to:

- Minimize inadvertent needlestick injuries with a passive safety needle-shielding feature.
- Improve catheter stability with an integrated stabilization platform designed to minimize movement within the vessel.
- Control the flow of blood which aids in the prevention of blood exposure.

The catheter may be used for any patient population with consideration given to adequacy of vascular anatomy and appropriateness for the solution being infused and duration of therapy.

Instructions for use

Materials used

Introcath Safety® 3 is available in two materials: PUR, ABS, Chrome-Nickel, Steel, POM, PP, Elastomer
FEP, ABS, Chrome-Nickel, Steel, POM, PP, Elastomer

Indications

- The Introcath Safety® 3 Closed Intravascular Catheter is inserted into a patient's vascular system for short term use (less than 30 days) to sample blood, monitor blood pressure, or administer fluids and blood intravascularly.
- The Introcath Safety® 3 is also indicated for subcutaneous infusion therapies.
- The 18-24 gauge catheters may be used with power injectors at a maximum pressure of 300 psi and a maximum flow rate. Use luer lock connection only.
See table below for recommended maximum flow rates for the Introcath Safety® IV Catheter when using a power injector. Testing performed at room temperature (22°C / 72°F). Warming contrast media per manufacturers' recommendations, can reduce the pressure necessary to achieve recommended flow rates.

	Contrast Media Viscosity [mPa's]	PUR
		Flow Rate (ml/sec)
24 Gauge*	2,3	5
	27,5	2,5
22 Gauge*	2,3	8
	27,5	5
20 Gauge*	2,3	14
	27,5	10
18 Gauge*	2,3	19
	27,5	15

* applicable to all lengths

Contraindications

Introcath Safety® 3 Closed IV Catheter should not be used in patients with known hypersensitivity to any of the materials employed.

Guidelines

- Observe Standard Precautions on ALL patients. Aseptic technique, proper skin preparation, and continued protection of the insertion site are essential.
- Use of only ISO Luer-slip and Luer-lock connectors ensures proper functioning.
- This device is intended for single use only and is provided sterile and non-pyrogenic.
- This device is sterile unless package has been opened or damaged.

- Components do not contain natural rubber latex, PVC, or DEHP.

Risks & Warnings

- Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.
- In the unlikely event that the safety mechanism does not engage, keep the needlepoint away from the body and fingers at all times and immediately dispose into an approved sharps container.
- Report needlesticks immediately and follow established protocols.
- In the case of an unsuccessful IV start, remove the needle first to activate safety mechanism, then remove catheter from patient and discard both.
- Do not bend the catheter/needle during insertion, advancement, or removal of the needle.
- Extreme care should be taken not to cut the catheter and possibly cause an embolism.
- Do not use scissors or sharp instruments at or near the insertion site.
- Never reinsert the needle inside the catheter once the needle has been partially or completely withdrawn as it may pierce and/or sever the catheter.
- Luer-slip connections should not be left unattended due to potential of disconnection.
- If blood is present, flush device according to institutional protocol.
- Do not attempt to override or defeat the safety shield.
- Patency of the catheter must be assured immediately before power injecting.
- Measures should be taken to avoid kinking or obstructing the catheter system during power injection to avoid product failure.

- Not intended to be used to facilitate the placement of Vascular Access devices such as guide-wires, indwelling central venous catheters, peripherally inserted central catheters and midline catheters into the vascular system.

Duration of use

- Change device according to CDC Guidelines and / or organizational policies.
- The insertion site should be checked at regular intervals.
- Introcan Safety® 3 Closed IV Catheter should be removed in the event of local or systemic signs of infection.

Application

1. Confirm all items are accessible for Introcan Safety® 3 Closed IV Catheter insertion/stabilization procedure.
2. Select and prepare site per organizational policies.
3. Apply tourniquet.
4. Holding flashback chamber, remove protective guard in a straight outward motion. Inspect device.
5. Do not rotate catheter hub prior to insertion.
6. Confirm catheter hub is properly seated on flashback chamber.
7. Hold the device by the flashback chamber with thumb and fingers on opposite sides. Verify push-off plate and needle bevel are in the "up" position.
8. Anchor the vessel with gentle skin traction. Adjust angle of insertion and access vessel. Observe blood return in flashback chamber to confirm vessel entry (see Fig. A).
9. Lower and advance the entire catheter and needle unit slightly to ensure catheter tip is in vessel (see Fig. B).
10. Using the push-off plate, advance catheter off the needle (approximately 1/8" or 3mm) and observe blood return between the needle

and catheter to confirm catheter is in the vessel (see Fig. C). After confirmation, continue advancing catheter off the needle and into vessel.

11. Release the tourniquet.
12. Gently press the stabilization platform to the skin to stabilize the catheter (see Fig. D).
13. Withdraw needle in a controlled and continuous motion. The safety shield automatically covers needle bevel as it exits the catheter hub (see Fig. E).
14. Immediately discard shielded needle into a puncture resistant, leak proof sharps container.
15. Immediately connect to tubing or adapter. Secure and apply sterile dressing (see Fig. F).

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B. Braun Medical AG
B. Braun Melsungen AG