
Discofix® C

Type Safeflow


en	Instructions for use	4
de	Gebrauchsanweisung	5
bg	Инструкция за употреба	7
cs	Návod k použití	9
da	Brugsanvisning	11
el	Οδηγίες χρήσης	13
es	Instrucciones de uso	15
et	Kasutusjuhised	16
fi	Käyttöohje	18
fr	Mode d'emploi	20
hr	Upute za upotrebu	22
hu	Használati utasítás	24
id	Petunjuk penggunaan	25
it	Istruzioni per l'uso	27
lt	Naudojimo instrukcija	29
lv	Lietošanas instrukcija	31
nl	Gebruikersinformatie	32
no	Bruksanvisning	34
pl	Instrukcja użytkowania	36
pt	Instruções de utilização	38
pt BR	Instruções de uso	40
ro	Instrucțiuni de utilizare	41
ru	Инструкция по применению	43
sk	Návod na použitie	45
sl	Navodila za uporabo	47
sr	Uputstvo za upotrebu	49
sv	Bruksanvisning	51

B | BRAUN

LLDorder 7558 - Grunewald 240735

Operating Instructions

en Instructions for use

Device description

Discofix® C stopcocks and manifolds with integrated valve (Discofix® C Type Safeflow) can be used with devices which are compatible or in conformance with Luer connectors, (e.g. ISO 80369-7). For gravity or pressure use up to 2 bar.

Sterility

Sterilized using gamma irradiation or ethylene oxide. Please refer to the label of primary packaging.

Intended purpose

Discofix® C stopcock Type Safeflow: Multidirectional stopcock for infusion therapy and monitoring.

Discofix® C stopcock Type Safeflow with extension line: Multidirectional stopcock with tubing for infusion therapy.

Discofix® C manifold Type Safeflow: Stopcock manifold for infusion therapy and monitoring.

Discofix® C manifold with extension line, Type Safeflow: Stopcock manifold with tubing for infusion therapy.

Indication

Discofix® C Type Safeflow is intended to be used for infusion and injection therapy, blood sampling and aspiration of body fluids as clinically indicated according to the Summary of product characteristics of the drugs / solutions. Variants without tube can also be used as control element for invasive pressure monitoring.

Patient population

Discofix® C Type Safeflow is intended to be used for all patients for whom infusion therapy is prescribed and / or for whom invasive pressure monitoring is necessary.

Intended User

Discofix® C Type Safeflow should only be used by authorized persons, e.g. healthcare professionals, licensed and / or certified physicians, nurses (according to local regulations) who have been adequately trained within the educational background in this technique. Upon medical assessment and adequate instruction, patients and/or caregivers may also be allowed to take over defined han-

dling steps according to national guidelines. Healthcare workers are responsible for ensuring that patients and / or caregivers are instructed in the correct use of Discofix® C Type Safeflow and that the instruction of caregiver and / or patient is documented by the healthcare worker.

Contraindications

Not known.

Residual Risks/Side effects

General risks occurring during infusion therapy include air embolism, local and systemic infection, particle contamination, medication error including over and under supply. Kinking of the line may result in decrease or stop of infusion flow. Leakage of drugs / fluids can be an additional hazard to patients, visitors or staff.

Precautions

- Visually inspect the sterile barrier system for breaches and integrity prior to use.
- Do not re-sterilize
- Flush product per national and / or institutional protocol.
- Mind the manufacturer's instructions for use of complementary used products (e.g. IV Set or pump).
- Please also consider the corresponding Summary of product characteristics of the drugs/ solutions regarding possible incompatibilities of drugs / solutions.

Warning

- Do not use if packaging or product is damaged or contaminated or if protective caps are loose or missing. Use of damaged or contaminated device might lead to injury, illness or death of the patient.
- Do not re-use: The products are single-use, sterile and non-pyrogenic medical devices. Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and / or impairment of functional capability. Contamination and / or limited functionality of the device may lead to injury, illness or death of the patient.
- Monitor all connections for tightness during the whole application.
- Do not re-attach the protective cap of the patient connector.

- Safeflow port: Luer slip syringes should be hand-tightened into the valve using a 90° turn. Do not leave luer slip syringes unattended.
- Safeflow port: Only use devices in conformance with Luer connectors (e.g. ISO 80369-7). Syringes or connectors which are not in conformance with Luer connectors (e.g. ISO 80369-7) can damage the swabable valve. Syringes and male luer connectors have a large variety of configurations and can vary significantly in design and dimensions.
- Do not use, in any situation, needles in the needle-free injection port.
- Stopcocks and manifolds should be primed completely. If a stopcock or a manifold is not primed carefully, up to 0.12 ml air per valve can remain in the infusion system (Total values: manifold 3-gang 0.36ml, manifold 5-gang 0.60ml).

Operating Instructions

Depending on the position of the handle the ports of the stopcock or manifold are open or closed. The handle knobs indicate that a port is open. A tactile response at every 45° turn confirms the exact positioning of the adjustment desired.

Mind institutional protocols and national guidelines for hygienic procedures such as disinfection and flushing during entire process.

Duration of use

Discofix® C Type Safeflow should be changed acc. to national standards (e.g. CDC) and / or institutional protocols. Mind that duration of use is dependent on the intended therapy according to the Summary of product characteristics of the solution or drug.

Disposal

Disposal acc. to local guidelines and/or clinical protocols.

Storage and handling conditions

Storage conditions:
Keep away from sunlight.
Keep dry.

Notice to the user

If, during the use of this product or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

If further Instructions for Use are required, they can be requested from the manufacturer or can be obtained from the B. Braun homepage: <https://eifu.bbraun.com/>.

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B. Braun Medical AG
B. Braun Melsungen AG