VICRYL RAPIDE™ Polyglactin 910 Suture

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Instructions for Use

VICRYL RAPIDE™ Polyglactin 910 Suture

Sterile Synthetic Absorbable Surgical Suture

Description

VICRYL RAPIDE™ Suture is a sterile, coated, braided, synthetic absorbable surgical suture composed of ≥ 87.8 weight percent Polyglactin 910, a copolymer made from 90 weight percent glycolide and 10 weight percent L-lactide. The empirical formula of the copolymer is $(C_2H_2O_2)_m(C_3H_4O_2)_n$. The characteristic rapid loss of strength is achieved by use of a polymer material with a lower molecular weight than regular VICRYL™ (Polyglactin 910) Suture.

VICRYL RAPIDE™ Suture has a coating with up to 12.0 weight percent of the suture composed of a copolymer made from 90 weight percent caprolactone and 10 weight percent glycolide and a mixture composed of approximately equal parts of calcium stearate and Polyglactin 370, a copolymer made from 30 weight percent glycolide and 70 weight percent lactide. Optionally, some codes are coated only with a mixture composed of approximately equal parts of Polyglactin 370 and calcium stearate. Polyglactin 910 copolymer and its coatings have been found to be nonantigenic, nonpyrogenic and elicit only a mild tissue reaction during absorption.

VICRYL RAPIDE™ Suture are available undyed as well as dyed with up to 0.20 weight percent D&C violet No. 2 (Colour Index Number 60725) to enhance visibility in the surgical field.

VICRYL RAPIDE™ Suture is available in a range of gauge sizes and lengths, non-needled or attached to stainless steel needles of varying types and sizes. Full details are listed in the catalogue.

VICRYL RAPIDE™ Suture complies with the requirements of the European Pharmacopoeia (Ph. Eur.) except for diameter for Sterile Synthetic Absorbable Braided Sutures and United States Pharmacopoeia (USP) for Absorbable Surgical Sutures, except for

a slight oversize in diameter for the gauges as shown in the table below and the exception of knot tensile strength. Knot tensile strength meets the USP and EP for collagen sutures (*Chorda resorbilis sterilis*).

USP Suture Size	Metric Suture Size	Maximum Oversize
6-0	0.7	0.008 mm
5-0	1	0.016 mm
4-0	1.5	0.017 mm
3-0	2	0.018 mm
2-0	3	0.010 mm
0	3.5	0.022 mm

The European Pharmacopoeia recognizes units of measure Metric and Ph. Eur. Sizes as equivalent, which is reflected on the labeling.

VICRYL RAPIDE™ Suture is intended for use only by healthcare professionals who are trained in surgical suturing techniques.

Clinical benefits expected from soft tissue approximation and where only short-term wound support is required are to promote wound healing by primary intention.

A summary of safety and clinical performance can be found at the following link (upon activation):

<https://ec.europa.eu/tools/eudamed>.

Indications / Intended Use

VICRYL RAPIDE™ Suture is intended for use in soft tissue approximation where only short-term wound support is required and where the rapid absorption of the suture would be beneficial. Due to its absorption profile VICRYL RAPIDE™ Suture is useful

for skin closure, particularly in paediatric surgery, episiotomies, circumcision, and closure of oral mucosa. VICRYL RAPIDE™ Suture is also successfully used in ophthalmic surgery for conjunctival sutures.

Contraindications

Due to the rapid loss of tensile strength, VICRYL RAPIDE™ Suture should not be used where extended approximation of tissues under stress is required or where wound support or ligation beyond 7 days is required. VICRYL RAPIDE™ Suture is not for use in cardiovascular and neurological tissues.

Warnings

Healthcare professionals should be familiar with surgical procedures and techniques involving absorbable sutures before employing VICRYL RAPIDE™ Suture for wound closure, as risk of tissue separation / wound dehiscence leading to impaired healing may vary with the site of application and the suture material used. Healthcare professionals should consider the *in vivo* performance (under the Performance / Actions section) when selecting a suture.

Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

As this is an absorbable suture material, the use of supplemental non-absorbable sutures should be considered by the healthcare professional in the closure of the sites that may undergo expansion, stretching or distension, or that may require additional support.

As an absorbable suture, this suture may act transiently as a foreign body. As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. Like all foreign bodies, VICRYL RAPIDE™ Suture may potentiate infection.

Do not resterilize/reuse. Reuse of this device (or portions of this device) may create a risk of product degradation that may result in device failure and/or cross-contamination that may lead to infection or transmission of blood-borne pathogens to patients and healthcare professionals.

Discard unintentionally opened / partially used / used devices and packages.

Precautions

As with any suture material, adequate knot security requires the standard surgical technique of flat, square ties with additional throws as warranted by surgical circumstance and the experience of the healthcare professional. Improper knotting or damage to suture during use can cause extended surgery time or treatment failure and additional surgery.

Skin sutures that must remain in place longer than 7 days may cause localized irritation and should be snipped off or removed as indicated by the healthcare professional. Subcuticular sutures should be placed as deeply as possible to minimize the erythema and induration normally associated with the absorption process.

Under some circumstances, notably orthopedic procedures, immobilization of joints by external support may be employed at the discretion of the healthcare professional.

Consideration should be taken in the use of absorbable sutures in tissues with poor blood supply as suture extrusion and delayed absorption may occur.

This suture may be inappropriate in elderly, malnourished or debilitated patients, or in patients suffering from conditions that may delay wound healing.

In handling this or any other suture material, care should be taken to avoid damage. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

Care should be taken to avoid damage when handling surgical needles. Grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the attachment end to the point. Grasping in the point area could impair the penetration performance and cause fracture of the needle. Grasping at the attachment end could cause bending or breakage. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking.

Healthcare professionals should exercise caution when handling surgical needles to avoid inadvertent needle-stick injury that may result in transmission of blood-borne pathogens from contaminated needles. Broken needles may result in extended or additional surgeries or residual foreign bodies. In the event of a product malfunction before use, such as a bent, broken or detached needle, or suture damage, the product should be discarded and a new one obtained to begin the procedure. In the event of product malfunction during use, it is up to the discretion of the healthcare professional whether to continue or discontinue usage of the product and how to complete the procedure.

Adverse Reactions / Undesirable Side Effects

Adverse reactions associated with the use of this device include tissue separation / wound dehiscence leading to impaired healing or failure to provide adequate wound support in the closure of the sites where expansion, stretching, or distension occur, and in elderly, malnourished or debilitated patients or in patients suffering from other conditions that may delay wound healing. Other associated adverse events include calculus formation in urinary or biliary tracts when prolonged contact with salt solutions such as urine or bile occurs, transient inflammatory tissue

reaction, and delayed absorption in tissue with poor blood supply. Transient local irritation and suture extrusion may occur at the wound site when skin sutures are left in place for greater than 7 days. Like all foreign bodies, VICRYL RAPIDE™ Suture may potentiate infection.

Suture breakage may result in hemorrhage. Suture material cutting or tearing through the tissue at the time of suture tensioning may cause soft tissue injury. Improper knotting or damage to suture during use can cause extended surgery time or treatment failure and additional surgery.

Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle-sticks with contaminated surgical needles may result in the transmission of blood-borne pathogens.

Healthcare professionals should convey adverse reactions, undesirable side effects and risks associated with the product and the procedure to the patient and advise the patient to contact a healthcare professional in case of any deviation from the normal postoperative course.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the country-competent authority.

Magnetic Resonance Imaging (MRI) / Carcinogenic, Mutagenic, Toxic to Reproduction (CMR) / Endocrine Disrupting (ED) Safety Information

VICRYL RAPIDE™ Suture is MR safe. No known CMR Category 1a/1b and ED substances are present at >0.1%. Category 1a/1b are defined as known or presumed human carcinogen (H340), mutagen (H350) or reproductive toxicant (H360), based on human evidence and animal studies.

Application / Instructions for Use

Sutures should be selected and implanted depending on patient condition, surgical experience, surgical technique, and wound characteristics. Discard needles in “sharps” containers. Discard unintentionally opened / partially used / used devices and packages according to your facility’s policies and procedures concerning biohazardous materials and waste.

VICRYL RAPIDE™ Suture typically falls off 7-10 days post-operatively or can be wiped off subsequently with sterile gauze. Normally the removal of the suture is not required.

Performance / Actions

VICRYL RAPIDE™ Suture elicits an initial, slight inflammatory reaction in tissue and ingrowth of fibrous connective tissue. Progressive loss of tensile strength and eventual absorption occurs by means of hydrolysis, where the copolymer degrades to

glycolic and lactic acids, which are subsequently absorbed and metabolized in the body. Absorption begins as a loss of tensile strength followed by a loss of mass.

All of the original tensile strength is lost by approximately 10 to 14 days post-implantation. Absorption is complete by 42 days. Implantation studies in rats indicate that VICRYL RAPIDE™ Suture retains the percentage of original tensile strength as follows:

Days Implantation	Approximate Original Strength Remaining
5 Days	50%

Sterility

VICRYL RAPIDE™ Suture that is dyed is sterilized by ethylene oxide gas and VICRYL RAPIDE™ Suture that is undyed is sterilized by irradiation. Do not resterilize. Do not use if package is opened or damaged.

Storage

No special storage required. Do not use after expiry date.

How Supplied

Please note that not all sizes are available in all markets. Please contact your local sales representative for size availability.

VICRYL RAPIDE™ Suture is available as sterile strands in sizes 8-0 through 1 (metric sizes 0.4-4.0), in a variety of lengths, with and without permanently attached needles.

VICRYL RAPIDE™ Suture is available in one, two or three dozen units per box.

Traceability

The following specific information can be found on the device packaging label: Catalogue number, Batch Code, expiry and manufacturing date, manufacturer name, address and the website and a Unique Device Identification bar code with the Global Trade Item Number information.

Symbols Used on Labeling

Catalogue number

Medical Device

Do not use if package is damaged and consult instructions for use

Caution

Do not re-use

Do not re-sterilize

Single sterile barrier system with protective packaging inside

Sterilized using ethylene oxide

Sterilized using irradiation

Batch code

Date of manufacture

Use-by date

Unique Device Identifier

Manufacturer

Authorized representative in the European Community/European Union

Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

Packaging unit

PEEL/LIFT/OPEN HERE

TEAR HERE

www.e-ifu.com
EU: +800 8888 2020
EU: +32 2 4037222

Consult instructions for use or consult electronic instructions for use.
EU: Call paper on demand help desk to get paper copies free of charge within 7 days.

MR Safe

Dyed - Absorbable - Braided - Coated - Suture

Undyed - Absorbable - Braided - Coated - Suture

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AB: Basılı kopyaları 7 gün içinde ücretsiz olarak talep üzerine almak için danışma hattını arayın.
Ethicon, LLC
Johnson & Johnson Medical GmbH