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### Class I medical devices

#### INSTRUCTIONS FOR PREPARATION

This procedure applies for all class I medical devices MEDIN, a.s. (hereinafter referred to as "devices"), except those directly intended for sterilization by the manufacturer before use. If the products are to be cleaned, disinfected and/or sterilized before use, follow the procedure below. The procedure for preparation a Class I device before use is at the discretion of the user.

#### 1. Caution

- 1.1. When cleaning and disinfecting medical devices made of stainless steel, it is necessary to avoid cleaning agents containing chloride ions, which can cause corrosion of steel parts.
- 1.2. Do not put the stainless steel products in the saline (NaCl), because the prolonged contact leads to the corrosion.
- 1.3. Hydrogen peroxide-based solutions and/or high-alkaline cleaners may cause discoloration of titanium parts.
- 1.4. Do not use abrasive cleaning materials such as scouring pad or abrasive detergents for cleaning. If possible, clean the devices in the washer/disinfector according to ISO 15883.

#### 2. Limits on reprocessing

- 2.1. No special precaution required.

#### 3. Treatment at the point of use

- 3.1. No special precaution required.

#### 4. Transportation

- 4.1. No special precaution required.

#### 5. Preparation for processing

- 5.1. No special precaution required.

#### 6. Automatic processing method

- 6.1. A validated automatic cleaning and disinfecting method must always be preferred over manual method for greater process safety.
- 6.2. It is recommended to wash devices before sterilisation using automatic washers with the option of thermal disinfection.
- 6.3. Place the devices in the washer so that they cannot move freely or overlap during the washing process and so that their surface cannot be damaged.
- 6.4. The use of a specific cleaning and disinfecting agent is not prescribed.
- 6.5. When using a washing and disinfecting agent, follow the product manufacturer's recommendations, observing the prescribed concentration, solution temperature, water quality and exposure time.

#### 7. Manual Processing method

- 7.1. Wash devices manually, under running water and then in a fresh cleaning and disinfecting solution in an ultrasonic bath, then rinse again under running water.
- 7.2. The use of a specific cleaning and disinfecting agent is not prescribed.
- 7.3. When using a washing and disinfecting agent, follow the product manufacturer's recommendations, observing the prescribed concentration, solution temperature, water quality and exposure time. High concentrations or prolonged exposure may cause damage to the device.
- 7.4. Remove the device from the cleaning solution after the recommended exposure time and rinse under running water.

#### 8. Disinfection

- 8.1. In case of cleaning in washer, disinfect at a minimum temperature of 90 °C for at least 5 minutes (according to ISO 15883-1).

#### 9. Drying

- 9.1. If the automatic washer does not include a drying cycle: dry each device thoroughly.

#### 10. Inspection and Maintenance

- 10.1. Before sterilization and before use, check that the device is undamaged and functional.

#### 11. Packaging

- 11.1. If the devices are to be sterilized before use, pack them separately in sterilization containers that meet the requirements of ISO 11607.

#### 12. Sterilization

- 12.1. The sterilization of the devices is not necessary.
- 12.2. If sterilization is required, sterilize the packaged devices with moist heat in a steam sterilizer at a temperature of 134 °C and a sterilization time of at least 5 minutes and a pressure of 304 kPa.
- 12.3. Only cleaned and disinfected devices may be sterilized.

#### 13. Storage

- 13.1. Store packaged devices in a dry, clean environment out of direct sunlight.

#### 14. Additional Information

- 14.1. The cleaning process with thermal disinfection has been verified by MEDIN, a.s. The manufacturer cannot guarantee the result if washing method other than those recommended is used. If a different cleaning method is used, it is recommended that this procedure be validated by the supplier of the particular washing equipment.
- 14.2. Instructions for cleaning and sterilization are provided according to ISO 15883.
- 14.3. If a suspected serious adverse event occurs in connection with the use of a medical device, report it to the manufacturer and to the national competent authority responsible for vigilance. You can find contact information for the manufacturer at [www.medin.cz](http://www.medin.cz) or communicate the suspicion directly to your contact person at MEDIN, a.s.

#### 15. Disposal

- 15.1. All packaging materials that are contaminated with human blood, tissues, secretions or body fluids and discarded products present a biological hazard and are considered as infectious waste. Dispose of this material in accordance with local regulations for medical infection.
- 15.2. Dispose of demonstrably uncontaminated material, such as uncontaminated packaging, in accordance with local regulations for medical waste.

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